




## Genesee Midcab Retractor



- Reversible rack with hinged arms and interchangeable, plug-in blades allow great versatility in retractor positioning.
- Adjustable angle blades incline rack in MI DCAB surgery.
- Deep blade retracts lung in minimally invasive mitral surgery.
- Excellent exposure in Mitral, LAD and IMA procedures.
- Unique ball clamps, rotary H-bar and pillar clamp allow fast setup and minimize snagging of sutures.
- Catalog number MIRS-1 (without attachments) and MIRS-AR (with attachments)



**Genesee BioMedical, Inc.**  
*Peak Performance Products for Cardiac Surgery*

1.800.786.4890

[www.geneseebiomedical.com](http://www.geneseebiomedical.com)