


genesee

## Intra-Art™ Aortic Punch

Reusable aortic punch with  
interchangable and  
reusable tips



- High Quality Stainless Steel Precision Instrument
- Reusable Punch with Easily Interchangeable Tips
- Extremely Cost Effective Compared to Disposable Punches



**Genesee BioMedical, Inc.**

*Peak Performance Products for Cardiac Surgery*

1.800.786.4890

[www.geneseebiomedical.com](http://www.geneseebiomedical.com)



## Intra-Art™ Multi-Size, Reusable Aortic Punch

- Easily Disassembled for Cleaning
- Coated Tips Inhibit Wear
- Each Tip Yields Approximately 250 or More Clean Precise Punches\*
- Custom Sterilization Tray Available for Punch and Tips
- Tips Available in Eight Different Sizes:
  - 3.0mm      3.6mm
  - 4.0mm      4.4mm
  - 4.8mm      5.2mm
  - 5.6mm      6.3mm
- Punch and Tips can be Purchased Individually or as a Complete Set



### Ordering Information

Catalog Number	Product Description
APCS	Aortic Punch with Complete Set of Tips and Sterilization Tray
APBT	Aortic Punch with One Tip (any size)
APT-3.0	3.0mm Aortic Punch Tip
APT-3.6	3.6mm Aortic Punch Tip
APT-4.0	4.0mm Aortic Punch Tip
APT-4.4	4.4mm Aortic Punch Tip
APT-4.8	4.8mm Aortic Punch Tip
APT-5.2	5.2mm Aortic Punch Tip
APT-5.6	5.6mm Aortic Punch Tip
APT-6.3	6.3mm Aortic Punch Tip
APST	Sterilization Tray

*\*Must follow instructions for use to yield greatest punches per tip.*

